

TILT EASE III CHEST - SPINE - HEAD

- ⊙ Solid Construction
- ⊙ Multi-Purpose
- ⊙ 130 KVP Power
- ⊙ Digital Displays
- ⊙ Touch Control Positioning
- ⊙ Electromagnetic And Disc Tilt Brakes
- ⊙ Precision Alignment
- ⊙ Adjustable Carriages
- ⊙ Powder Coat Paint
- ⊙ Cost Effective

⊙ Rotates 360°

⊙ Angles -45° to +90°



TXR
SYSTEM 4

SYSTEM FOUR

TILT EASE III – CHEST-SPINE-HEAD FULL TILT SYSTEM "TOP OF THE LINE"

Equipped with 14" x 17" GRIDCABINET

Product No.	Description	Suggested List
325 D	Standard Frequency 300 mA, 130 kVp X-Ray Control with kVp, mA, & mAs digital readouts... <i>Provided with SCR's, tube limit protection and self-diagnostic program.</i>	\$ 15,625.00
X-Ray Tube	Toshiba 1.0-2.0 mm focus tube E7239FX 90° or Toshiba 1.0-2.0 mm focus tube E7239GX 180°	
13-SM-84	Tubestand with magnetic locks, 6' floor rail and Analog Positioning Control	
24-15500	Proflex 15' High Tension cables (<i>pair</i>)	
70-63000-1	Certified Manual Collimator	
GC-NS	14" x 17" Gridcabinet with Heavy Duty Tray , Panel Front and 8:1 103 line grid . Focus 40-72 inches.	
13-BT-7 III	Heavy Duty Precision Manual Tilt Ease. <i>(Accepts Buckys or Gridcabinets. Angle scale and lock are manual, vertical lock is electromagnetic)</i>	

OPTIONS

High Frequency	High Frequency Options: Deduct the 325D List Price from the System above, then Add the High Frequency List Price.....	For Pricing See Page 9
HC	Head Clamp	\$ 990.00
2341-196	Chair, Swivel – Rotary with floor pad (See System Five).....	300.00
E7242FX	Toshiba 90° .6 - 1.5 mm fractional focus tube	1,170.00
E7242GX	Toshiba 180° .6 - 1.5 mm fractional focus tube	1,170.00
325M	Standard Frequency 300 mA X-Ray Control.....	300.00
525 SFQ	Standard Frequency 500 mA X-Ray Control.....	900.00
50Hz Option	50 Hertz Stabilizer – Includes Timer Dial.....	200.00
24-20500	Proflex 20' High Tension cables (<i>pair</i>).....	75.00
17ANS	14" x 17" Par Speed Bucky 8:1 103 line grid (<i>each</i>)	1,125.00
10-08010	14" x 17" Super Speed Bucky 8:1 103 line grid (<i>each</i>)	1,370.00
1011034072	Upgrade grid to 18 x 18 10:1 103 line 40-72" (<i>each</i>)	50.00

To add Automatic Exposure Control to Single phase or to order Ion Chambers, see page 3 in this booklet for pricing.

NOTES

Please check our web site for the latest updates to our product line at www.txr.com

Sample Room Layout for TXR System Four Shown with 325D X-Ray Control

